

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse IgG2a FITC Negative Control NMB-02A2

Artikelname	Mouse IgG2a FITC Negative Control
Artikelnummer	NMB-02A2
Hersteller Artikelnummer	02A2
Alternativnummer	NMB-02A2
Hersteller	NordicMubio
Wirt	Mouse
Kategorie	Antikörper
Applikation	FC
Konjugation	FITC
Produktbeschreibung	Negative control for immunofluorescence staining with FITC conjugated mouse monoclonal antibodies. Assessment of non-specific binding of mouse monoclonal antibodies to human cell surface antigens....
Konzentration	See vial for concentration
Isotyp	IgG2a
Puffer	Purified mouse IgG2a conjugated to FITC, as solution in phosphate buffered saline with 0.08% sodium azide and 0.2% carrier protein
Reinheit	Protein A/G Chromatography
Formel	Provided as solution in phosphate buffered saline with 0.08% sodium azide and 0.2% carrier protein

Anwendungsbeschreibung	Consult the appropriate fact sheet to determine the amount of antibodies to be used as a control for PBMC or Whole Blood. PBMC: Add 10 µl of MAB/10⁶ PBMC in 100 µl PBS. Mix gently and incubate for 15 minutes at 2 to 8°C. Wash twice with PBS and analyze. WHOLE BLOOD: Add 10 µl of MAB/100 µl of whole blood. Mix gently and incubate for 15 minutes at room temperature 20°C. Lyse the whole blood. Wash once with PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.
------------------------	--