

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

LPS from Salmonella abortus equi S-form (TLRGRADE) (Ready-to-Use) (Biotin) **ENZ-ALX-581-150-R500**

Artikelname	LPS from Salmonella abortus equi S-form (TLRGRADE) (Ready-to-Use) (Biotin)
Artikelnummer	ENZ-ALX-581-150-R500
Hersteller Artikelnummer	ALX-581-150-R500
Alternativnummer	ENZ-ALX-581-150-R500-500
Hersteller	Enzo Life Sciences
Kategorie	Sonstiges
Konjugation	Biotin
Produktbeschreibung	Useful reagent for evaluating LPS binding and cell activation in white blood cells....
Quelle	Smooth (S)-form LPS, isolated and purified from Salmonella abortus equi by modification of the phenol water extraction and PCP method. LPS-Biotin was prepared using the biotin reagent biotinamidocaproate N-hydroxysuccinimide ester. Briefly, LPS at 10mg/ml
Reinheit	Absence of detectable protein or DNA contaminants with agonistic TLR activity.
Formulierung	Liquid. Sterile, ready-to-use solution in pyrogen-free double-distilled water.
Target-Kategorie	TLR4