

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Y-79 (human retinoblastoma) nuclear extract lysate (denatured), Human ABN-L042V3

Artikelname	Y-79 (human retinoblastoma) nuclear extract lysate (denatured), Human
Artikelnummer	ABN-L042V3
Hersteller Artikelnummer	L042V3
Alternativnummer	ABN-L042V3-50
Hersteller	Abnova
Wirt	Human
Kategorie	Zellen/Zellkultur
Applikation	WB
Spezies Reaktivität	Human
Produktbeschreibung	Nuclear extract cell lysate (denatured)....
Puffer	In ready-to-use 1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue).
Quelle	Retina
Application Verdünnung	The lysate is ready to load on SDS-PAGE for Western blotting.
Anwendungsbeschreibung	Nuclear extract was prepared by using a modified protocol of Dignam et al. Cells were Harvested and homogenized in Buffer A, and then centrifugated at 25,000 g for 20 minutes to remove cytoplasm and pellet the nuclei. The pellet was re-suspended in Buffe