

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Chlamydia trachomatis Virulence plasmid protein pGP3-D protein BYT-ORB704784

Artikelname	Chlamydia trachomatis Virulence plasmid protein pGP3-D protein
Artikelnummer	BYT-ORB704784
Hersteller Artikelnummer	orb704784
Alternativnummer	BYT-ORB704784-1,BYT-ORB704784-100,BYT-ORB704784-20
Hersteller	Biorbyt
Kategorie	Molekularbiologie
Spezies Reaktivität	Bacteria
Produktbeschreibung	Recombinant Chlamydia trachomatis Virulence plasmid protein pGP3-D...
Molekulargewicht	30.3 kDa
UniProt	P0CE18
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Chlamydia trachomatis
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	GNSGFYLYNTQNCVFADNIKVGQMTEPLKDQQIILGTTSTPVAAKMTASDGIS LTVSNNPSTNASITIGLDAEKAYQLILEKLGDQILGGIADTIVDSTVQDILDKITT DPSLGLLKAFNNFPITNKIQCNGLFTPRNIETLLGGTEIGKFTVTPKSSGSMFLV SADIIASRMEGGVVLALVREGDSKPYAISYGYSSGVPNLCSLRTRIINTGLTPTT YSLRVGGLESGVWVWNALSNGNDILGITNTSNVSF
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration