

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

HEK293 Whole Cell Lysate Insulin Stimulated, Unconjugated BYT-ORB348680

Artikelname	HEK293 Whole Cell Lysate Insulin Stimulated, Unconjugated
Artikelnummer	BYT-ORB348680
Hersteller Artikelnummer	orb348680
Alternativnummer	BYT-ORB348680-500
Hersteller	Biorbyt
Kategorie	Sonstiges
Applikation	SDS-PAGE, WB
Konjugation	Unconjugated
Produktbeschreibung	HEK293 Whole Cell Lysate Insulin Stimulated...
Konzentration	1.0 mg/ml
Puffer	Preservative: None. Stabilizer: 10% (v/v) Glycerol. 1X SDS-PAGE Sample Buffer (62.5 mM Tris HCl, 2% SDS, 10% Glycerol and 0.005% bromophenol blue, pH 6.8)
Quelle	Human
Reinheit	The cells were grown in Dulbeccos medium supplemented with 10% fetal bovine serum. Cells were washed with PBS and then incubated on ice in modified RIPA buffer, containing 150 mM sodium chloride, 50 mM Tris HCl, pH 7.4, 1 mM EDTA, 1.0% NP-40, 0.5% sodiu
Formulierung	Liquid (sterile filtered)
Application Verdünnung	WB: User Optimized

Anwendungsbeschreibung	ready-to-use lysates are especially prepared as positive controls for separation by SDS-PAGE and subsequent western blot analysis. Lysates are prepared in denaturing buffer WITHOUT dissociating agents (i.e. no 2-mercaptoethanol or dithiothreitol has been
------------------------	--