

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

HeLa Cell Nuclear Extract Nocodazole Stimulated, Unconjugated BYT-ORB348686

Artikelname	HeLa Cell Nuclear Extract Nocodazole Stimulated, Unconjugated
Artikelnummer	BYT-ORB348686
Hersteller Artikelnummer	orb348686
Alternativnummer	BYT-ORB348686-200
Hersteller	Biorbyt
Kategorie	Sonstiges
Applikation	SDS-PAGE, WB
Konjugation	Unconjugated
Produktbeschreibung	HeLa Cell Nuclear Extract Nocodazole Stimulated...
Konzentration	1.0 mg/ml
Puffer	Preservative: None. Stabilizer: 10% (v/v) Glycerol. 1X SDS-PAGE Sample Buffer (62.5 mM Tris HCl, 2% SDS, 10% Glycerol and 0.005% bromophenol blue, pH 6.8)
Quelle	Human
Reinheit	The cells were grown in DMEM supplemented with 10% FBS (Fetal Bovine Serum). Cells were treated with 0.2 µg/ml Nocodazole for 30 min. The lysate was prepared by first washing the cells in PBS. Washed cells were then incubated on ice in lysis buffer conta
Formulierung	Liquid (sterile filtered)
Application Verdünnung	WB: User Optimized

Anwendungsbeschreibung	ready-to-use nuclear extracts are especially prepared as positive controls for separation by SDS-PAGE and subsequent western blot analysis. Nuclear extracts are supplied in denaturing buffer without dissociating agents. Heat nuclear extract to 95C for 5
------------------------	---