

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

RH HB-EGF **SHD-100-180-50UG**

Artikelname	RH HB-EGF
Artikelnummer	SHD-100-180-50UG
Hersteller Artikelnummer	100-180-50UG
Alternativnummer	SHD-100-180-50UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Sonstiges
Applikation	Bioassay, SDS-PAGE
Spezies Reaktivität	Human
Produktbeschreibung	Heparin-binding EGF-like growth factor (HB-EGF) is a member of the epidermal growth factor (EGF) family and is expressed by monocytes and macrophages. HB-EGF is the predominant growth factor involved in epithelialization during wound healing. HB-EGF ...
Molekulargewicht	9.9
UniProt	Q99075
Puffer	10 mM Sodium Phosphate, pH 7.5
Quelle	E.coli
Reinheit	95
Formulierung	Lyophilized from a 0.2 µm filtered solution

Sequenz	MDLQEADLDL LRVTLSSKPQ ALATPNKEEH GKRKKKGKGL GKKRDPCLRK YKDFCIHGEC KYVKELRAPC CICHPGYHGERCHGLSL
Anwendungsbeschreibung	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.