

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

RH ENA 78 **SHD-100-185-20UG**

Artikelname	RH ENA 78
Artikelnummer	SHD-100-185-20UG
Hersteller Artikelnummer	100-185-20UG
Alternativnummer	SHD-100-185-20UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Sonstiges
Applikation	Bioassay, SDS-PAGE
Spezies Reaktivität	Human
Produktbeschreibung	Epithelial-derived neutrophil-activating peptide 78 (ENA 78), also known as CXCL5, is a chemokine that recruits neutrophils, promotes angiogenesis, and stimulates connective tissue remodelling. ENA 78 production is stimulated by interleukin 1 (IL-1) ...
Molekulargewicht	8.1
UniProt	P42830
Puffer	0.1% Trifluoroacetic Acid (TFA)
Quelle	E.coli
Reinheit	95
Formulierung	Lyophilized from a 0.2 µm filtered solution

Sequenz	AAVLRELRCV CLQTTQGVHP KMISNLQVFA IGPQCSKVEV VASLKNGKEI CLDPEAPFLK KVIQKILDGG NKEN
Anwendungsbeschreibung	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.