

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

RH SDF-1 A / CXCL12 SHD-100-20-500UG

Artikelname	RH SDF-1 A / CXCL12
Artikelnummer	SHD-100-20-500UG
Hersteller Artikelnummer	100-20-500UG
Alternativnummer	SHD-100-20-500UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Sonstiges
Applikation	Bioassay, SDS-PAGE
Spezies Reaktivität	Human
Produktbeschreibung	Stromal cell-derived factor-1 alpha (SDF-1 alpha), also called CXCL12a, is one of two SDF-1 splice variants made by a wide variety of cells upon stimulation by inflammatory cytokines such as TNF, IL-1, and LPS. SDF-1alpha signals through the G protei...
Molekulargewicht	8
UniProt	P48061
Puffer	0.1% Trifluoroacetic Acid (TFA)
Quelle	E.coli
Reinheit	95
Formulierung	Lyophilized from a 0.2 µm filtered solution

Sequenz	KPVSLSYRCP CRFFESHVAR ANVKHLKILN TPNCALQIVA RLKNNNRQVC IDPKLKWIQE YLEKALNK
Anwendungsbeschreibung	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.