

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

RH FLT-3 LIGAND (AF) SHD-100-21AF-10UG

Artikelname	RH FLT-3 LIGAND (AF)
Artikelnummer	SHD-100-21AF-10UG
Hersteller Artikelnummer	100-21AF-10UG
Alternativnummer	SHD-100-21AF-10UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Sonstiges
Applikation	Bioassay, SDS-PAGE
Spezies Reaktivität	Human
Produktbeschreibung	Fms-related tyrosine kinase 3 ligand (FLT-3 ligand) is a growth factor that regulates hematopoietic cell proliferation. FLT-3 ligand signalling is transmitted through the fms-related tyrosine kinase 3 (FLT-3) receptor. FLT-3 ligand promotes the long...
Molekulargewicht	17.7
UniProt	P49771
Puffer	10 mM Sodium Phosphate, 50 mM Sodium Chloride, pH 7.5
Quelle	E.coli
Reinheit	97
Formulierung	Lyophilized from a 0.2 µm filtered solution

Sequenz	MTQDCSFQHS PISSDFAVKI RELSDYLLQD YPVTVASNLQ DEELCGGLWR LVLAQRWMER LKTVAGSKMQ GLLERVNTEI HFVTKCAFQP PPSCLRFVQT NISRLQETS EQLVALKPWI TRQNFSRCLE LQCQPDSSL PPPWSPRPLE ATAPT
Anwendungsbeschreibung	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.