

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

RM FLT-3 LIGAND SHD-200-06-100UG

Artikelname	RM FLT-3 LIGAND
Artikelnummer	SHD-200-06-100UG
Hersteller Artikelnummer	200-06-100UG
Alternativnummer	SHD-200-06-100UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Sonstiges
Applikation	Bioassay, SDS-PAGE
Spezies Reaktivität	Mouse
Produktbeschreibung	Fms-related tyrosine kinase 3 ligand (FLT-3 ligand) is a growth factor that regulates hematopoietic cell proliferation. FLT-3 ligand signaling is transmitted through the fms-related tyrosine kinase 3 (FLT-3) receptor. FLT-3 ligand promotes the long-t...
Molekulargewicht	18.6
UniProt	P49772
Puffer	10 mM Sodium Phosphate, pH 7.5
Quelle	E.coli
Reinheit	95
Formulierung	Lyophilized from a 0.2 µm filtered solution

Sequenz	MTPDCYFSHS PISSNFKVKF RELTDHLLKD YPVTVAVNLQ DEKHCKALWS LFLAQRWIEQ LKTVAGSKMQ TLLLEDVNTEI HFVTSCTFQP LPECLRFVQT NISHLLKDTC TQLLALKPCI GKACQNFSRC LEVQCQPDSS TLLPPRSPIA LEATELPEPR PRQ
Anwendungsbeschreibung	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.