

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### RR EGF SHD-300-05-20UG

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | RR EGF                                                                                                                                                                                                                                                        |
| Artikelnummer            | SHD-300-05-20UG                                                                                                                                                                                                                                               |
| Hersteller Artikelnummer | 300-05-20UG                                                                                                                                                                                                                                                   |
| Alternativnummer         | SHD-300-05-20UG                                                                                                                                                                                                                                               |
| Hersteller               | Fujifilm Irvine Scientific                                                                                                                                                                                                                                    |
| Kategorie                | Sonstiges                                                                                                                                                                                                                                                     |
| Applikation              | Bioassay, SDS-PAGE                                                                                                                                                                                                                                            |
| Spezies Reaktivität      | Rat                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | Epidermal growth factor (EGF) is a growth factor that stimulates the proliferation, differentiation, and survival of epithelial and epidermal cells. EGF contains three intramolecular disulfide bonds and binds in high affinity to the epidermal growth ... |
| Molekulargewicht         | 6.3                                                                                                                                                                                                                                                           |
| UniProt                  | <a href="#">P07522</a>                                                                                                                                                                                                                                        |
| Puffer                   | 10 mM Sodium Phosphate, pH 7.5                                                                                                                                                                                                                                |
| Quelle                   | E.coli                                                                                                                                                                                                                                                        |
| Reinheit                 | 95                                                                                                                                                                                                                                                            |
| Formulierung             | Lyophilized from a 0.2 µm filtered solution                                                                                                                                                                                                                   |

|                        |                                                                                                                                                                                                                                                                |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sequenz                | MNSNTGCPPS YDGYCLNGGV CMYVESVDY VCNCVIGYIG<br>ERCQHRDLRW WKLR                                                                                                                                                                                                  |
| Anwendungsbeschreibung | Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0. |