

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### RR IL-17F SHD-300-28-1MG

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | RR IL-17F                                                                                                                                                                                                                                                     |
| Artikelnummer            | SHD-300-28-1MG                                                                                                                                                                                                                                                |
| Hersteller Artikelnummer | 300-28-1MG                                                                                                                                                                                                                                                    |
| Alternativnummer         | SHD-300-28-1MG                                                                                                                                                                                                                                                |
| Hersteller               | Fujifilm Irvine Scientific                                                                                                                                                                                                                                    |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                              |
| Applikation              | Bioassay, SDS-PAGE                                                                                                                                                                                                                                            |
| Spezies Reaktivität      | Rat                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | Interleukin 17F (IL-17F), a member of the IL-17 cytokine family, is secreted by activated CD4+ T cells and monocytes. IL-17F binds the IL-17 receptor related molecule, IL17RC, to promote the production of the interleukin 6 (IL-6), interleukin 8 (IL-8... |
| Molekulargewicht         | 15.0/30.1                                                                                                                                                                                                                                                     |
| UniProt                  | <a href="#">Q5BJ95</a>                                                                                                                                                                                                                                        |
| Puffer                   | 0.1% Trifluoroacetic Acid (TFA)                                                                                                                                                                                                                               |
| Quelle                   | E.coli                                                                                                                                                                                                                                                        |
| Reinheit                 | 95                                                                                                                                                                                                                                                            |
| Formulierung             | Lyophilized from a 0.2 µm filtered solution                                                                                                                                                                                                                   |

|                        |                                                                                                                                                                                                                                                                            |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sequenz                | MARRNPKVGL SALQKAGNCP PLEDNSVRVD IRIFNQNQGI SVPRDFQNRS<br>SSPWDYNITR DPDRFPSEIA EAQCRHSGCI NAQGQEDGSM NSVPIQQEIL<br>VLRREPQGCS NSFRLEKMLI KVGCTCVTPI VHHA                                                                                                                  |
| Anwendungsbeschreibung | Centrifuge vial before opening. Suspend the product by gently<br>pipetting the above recommended solution down the sides of the<br>vial. DO NOT VORTEX. Allow several minutes for complete<br>reconstitution. For prolonged storage, dilute to working aliquots in a<br>0. |