

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

RR IL-17AF SHD-300-29-500UG

Artikelname	RR IL-17AF
Artikelnummer	SHD-300-29-500UG
Hersteller Artikelnummer	300-29-500UG
Alternativnummer	SHD-300-29-500UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Proteine/Peptide
Applikation	Bioassay, SDS-PAGE
Spezies Reaktivität	Rat
Produktbeschreibung	Interleukin 17AF (IL-17AF) is a heterodimer that is composed of the interleukin 17A (IL-17A) and interleukin 17F (IL-17F) members of the IL-17 family of cytokines. IL-17AF is produced by T helper 17 cells (Th17) following interleukin 23 (IL-23) stimu...
Molekulargewicht	30.2
UniProt	Q61453
Puffer	0.1% Trifluoroacetic Acid (TFA)
Quelle	E.coli
Reinheit	95
Formulierung	Lyophilized from a 0.2 µm filtered solution

Sequenz	IL-17A: MAVLIPQSSV CPNAEANNFL QNVKVNKVL NSLSSKASSR RPSDYLNIRST SPWTLNRNED PDRYPSVIWE AQCRHQRCVN AEGKLDHHMN SVLIQQEILV LKREPEKCPF TFRVEKMLVG VGCTCVSSIV RHAS IL-17F:MARRNPKVGL SALQKAGNCP PLEDNSVRVD IRIFNQNQGI SVPRDFQNRS SSPWD YNITR DPDRFPSEIA EAQCRHSGCI NA
Anwendungsbeschreibung	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.