

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

EIF3E Peptide PRS-42-702P

Artikelname	EIF3E Peptide
Artikelnummer	PRS-42-702P
Hersteller Artikelnummer	42-702P
Alternativnummer	PRS-42-702P-0.1
Hersteller	ProSci
Kategorie	Proteine/Peptide
Applikation	Blocking
NCBI	3646
UniProt	P60228
Puffer	EIF3E peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL.
Formulierung	Lyophilized

Anwendungsbeschreibung	EIF3E Peptide is used for blocking the activity of EIF3E antibody (catalog number 42-702). The EIF3E antibody ships at 0.5 mg/mL and the EIF3E Peptide ships as a 100 µg pellet. When peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL, identical to that of the antibody. The best starting ratio would be 1:1 (which means molar excess of peptide relative to antibody when identical volumes are mixed). Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been pre-absorbed with the provided peptide.
------------------------	---