

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Topoisomerase I, Mitochondrial (TOP1MT) (TOP1MT/2883R), CF640R conjugate, 0.1mg/mL, Clone: [TOP1MT/2883R], Rabbit, Monoclonal BOT-BNC402883-100**

Artikelname	Topoisomerase I, Mitochondrial (TOP1MT) (TOP1MT/2883R), CF640R conjugate, 0.1mg/mL, Clone: [TOP1MT/2883R], Rabbit, Monoclonal
Artikelnummer	BOT-BNC402883-100
Hersteller Artikelnummer	BNC402883-100
Alternativnummer	BOT-BNC402883-100-100UL
Hersteller	Biotium
Wirt	Rabbit
Kategorie	Antikörper
Applikation	IHC
Spezies Reaktivität	Human
Immunogen	Recombinant full-length human TOP1MT protein
Konjugation	CF640R
Produktbeschreibung	DNA topoisomerases are nuclear enzymes that regulate the topological structure of DNA in eukaryotic cells by transiently breaking and rejoining DNA strands. Due to their roles in DNA replication, recombination, and transcription, DNA topoisomerases h...
Klonalität	Monoclonal
Konzentration	0.1 mg/mL

Klon-Bezeichnung	[TOP1MT/2883R]
Molekulargewicht	70 kDa
UniProt	<a href="#">Q969P6</a>
Puffer	PBS, 0.1% BSA, 0.05% azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (formalin): 0.5-1 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user