

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Eosinophil Peroxidase(EPO104), CF647 conjugate, 0.1mg/mL, Clone: [EPO104], Mouse, Monoclonal BOT-BNC470104-100**

|                          |   |
|--------------------------|---|
| Artikelname              | Eosinophil Peroxidase(EPO104), CF647 conjugate, 0.1mg/mL, Clone: [EPO104], Mouse, Monoclonal  |
| Artikelnummer            | BOT-BNC470104-100   |
| Hersteller Artikelnummer | BNC470104-100   |
| Alternativnummer         | BOT-BNC470104-100-100UL   |
| Hersteller               | Biotium   |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Spezies Reaktivität      | Human   |
| Immunogen                | Human eosinophils from a patient with hypereosinophilic syndrome  |
| Konjugation              | CF647   |
| Produktbeschreibung      | Peripheral blood granulocytes are classified into neutrophils, basophils and eosinophils according to the staining characteristics of their cytoplasmic granules. Granule proteins are released by physiologic and pharmacologic stimuli and play importan... |
| Klonalität               | Monoclonal  |
| Konzentration            | 0.1 mg/mL   |
| Klon-Bezeichnung         | [EPO104]  |
| Molekulargewicht         | 70 kDa (55 kDa-Heavy Chain, 15 kDa-Light Chain)   |

|                        |  |
|------------------------|--|
| UniProt                | <a href="#">P11678</a>   |
| Puffer                 | PBS, 0.1% BSA, 0.05% azide   |
| Quelle                 | Animal   |
| Anwendungsbeschreibung | For coating for ELISA, order Ab without BSA Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Optimal dilution and staining procedure for a specific application should be determined by user Recommended starting concentrations for titration are 1-2 ug/mL for most applications, or 1 ug/million cells/100 uL for flow cytometry |