

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

**gp100 / Melanosome / PMEL17 / SILV (Melanoma Marker)(HMB45 + PMEL/783), Biotin conjugate, 0.1mg/mL, Clone: [HMB45 PMEL/783], Mouse, Monoclonal  
BOT-BNCB0951-100**

|                          |   |
|--------------------------|---|
| Artikelname              | gp100 / Melanosome / PMEL17 / SILV (Melanoma Marker)(HMB45 + PMEL/783), Biotin conjugate, 0.1mg/mL, Clone: [HMB45 PMEL/783], Mouse, Monoclonal  |
| Artikelnummer            | BOT-BNCB0951-100  |
| Hersteller Artikelnummer | BNCB0951-100  |
| Alternativnummer         | BOT-BNCB0951-100-100UL  |
| Hersteller               | Biotium   |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Applikation              | IHC, WB   |
| Spezies Reaktivität      | Human   |
| Immunogen                | Extract of pigmented melanoma metastases from lymph nodes (HMB45), Recombinant human SILV protein (PMEL/783)  |
| Konjugation              | Biotin  |
| Produktbeschreibung      | By immunohistochemistry, this antibody specifically recognizes a protein in melanocytes and melanomas. This MAb reacts with junctional and blue nevus cells and variably with fetal and neonatal melanocytes. Intradermal nevi, normal adult melanocytes, ... |
| Klonalität               | Monoclonal  |

|                        |  |
|------------------------|--|
| Konzentration          | 0.1 mg/mL  |
| Klon-Bezeichnung       | [HMB45 PMEL/783]   |
| Molekulargewicht       | 90-100 kDa   |
| UniProt                | <a href="#">P40967</a>   |
| Puffer                 | PBS, 0.1% BSA, 0.05% azide   |
| Quelle                 | Animal   |
| Anwendungsbeschreibung | Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistochemistry (formalin-fixed): 0.5-1.0 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user |