

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) (CALD1/1424R), Biotin conjugate, 0.1mg/mL, Clone: [CALD1/1424R], Rabbit, Monoclonal BOT-BNCB1424-100**

Artikelname	Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) (CALD1/1424R), Biotin conjugate, 0.1mg/mL, Clone: [CALD1/1424R], Rabbit, Monoclonal
Artikelnummer	BOT-BNCB1424-100
Hersteller Artikelnummer	BNCB1424-100
Alternativnummer	BOT-BNCB1424-100-100UL
Hersteller	Biotium
Wirt	Rabbit
Kategorie	Antikörper
Applikation	IHC, WB
Spezies Reaktivität	Human, Rat
Immunogen	Recombinant human full-length CALD1 protein
Konjugation	Biotin
Produktbeschreibung	Caldesmon HMW is the high molecular weight variant of Caldesmon. Two closely related variants of human caldesmon have been identified which are different in their electrophoretic mobility and cellular distribution. The h-caldesmon variant (120-150 kD...
Klonalität	Monoclonal
Konzentration	0.1 mg/mL

Klon-Bezeichnung	[CALD1/1424R]
Molekulargewicht	150 kDa
UniProt	<a href="#">Q05682</a>
Puffer	PBS, 0.1% BSA, 0.05% azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunofluorescence: 1-2 ug/mL Immunohistology (formalin) 1-2 ug/mL Staining of formalin-fixed tissues requires boiling tissue sections in 1 mM EDTA pH 7.5-8.5 Flow Cytometry 0.5-1 ug/million cells/0.1 mL Optimal dilution for a specific application should be determined by user