

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Melanoma Marker(HMB45 + A103 + T311), 0.2mg/mL, Clone: [HMB45 A103 T311], Mouse, Monoclonal BOT-BNUB0702-100

Artikelname	Melanoma Marker(HMB45 + A103 + T311), 0.2mg/mL, Clone: [HMB45 A103 T311], Mouse, Monoclonal
Artikelnummer	BOT-BNUB0702-100
Hersteller Artikelnummer	BNUB0702-100
Alternativnummer	BOT-BNUB0702-100-100UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	IHC, WB
Spezies Reaktivität	Human
Immunogen	Recombinant hMART-1 protein (A103), Recombinant tyrosinase protein (T311), Extract of pigmented melanoma metastases from lymph nodes (HMB45)
Produktbeschreibung	This antibody cocktail recognizes three melanoma-specific proteins, which include MART-1, Tyrosinase and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Tyrosinase is one of the...
Klonalität	Monoclonal
Konzentration	0.2 mg/mL
Klon-Bezeichnung	[HMB45 A103 T311]

Molekulargewicht	20-22 kDa (doublet) (MART), 70-80 kDa (Tyrosinase), 90-100 kDa (gp100)
UniProt	Q16655
Puffer	PBS, 0.05% BSA, 0.05% azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (Formalin-fixed) 0.5-1.0 ug/mL Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user