

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat Anti-Rabbit IgG (H&L) Highly Cross-Adsorbed - AF488 BWT-BS22352

Artikelname	Goat Anti-Rabbit IgG (H&L) Highly Cross-Adsorbed - AF488
Artikelnummer	BWT-BS22352
Hersteller Artikelnummer	BS22352
Alternativnummer	BWT-BS22352-100UL,BWT-BS22352-500UL
Hersteller	Bioworld Technology
Wirt	Goat
Kategorie	Antikörper
Applikation	ELISA, FC, IF, IH
Immunogen	Rabbit IgG
Konjugation	AF488
Produktbeschreibung	Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) Highly Cross-Adsorbed and conjugate to AF488...
Konzentration	1 mg/ml
Puffer	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Reinheit	<p>Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) have been cross-adsorbed against IgG from bovine, goat, guinea pig, hamster, horse, mouse, rat, and human. Cross-adsorption or pre-adsorption is a purification step to increase specificity of the antibody resulting in less background staining and cross-reactivity. The secondary antibody solution is passed through a column matrix containing immobilized serum proteins from potentially cross-reactive species. Only the nonspecific-binding secondary antibodies are captured in the column, and the highly specific secondaries flow through. Further passages through additional columns result in highly cross-adsorbed preparations of secondary antibody. The benefits of these extra steps are apparent in multiplexing/multicolor-staining experiments where there is potential cross-reactivity with other primary antibodies or in tissue/cell fluorescent staining experiments where there may be the presence of endogenous immunoglobulins.</p>
Application Verdünnung	<p>IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)</p>