

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

CTGRADE GMP RH M-CSF SHD-500-12-50UG

Artikelname	CTGRADE GMP RH M-CSF
Artikelnummer	SHD-500-12-50UG
Hersteller Artikelnummer	500-12-50UG
Alternativnummer	SHD-500-12-50UG
Hersteller	Fujifilm Irvine Scientific
Kategorie	Sonstiges
Applikation	Bioassay, SDS-PAGE, WB
Spezies Reaktivität	Human
Produktbeschreibung	Macrophage colony-stimulating factor (M-CSF) is a hematopoietic growth factor with broad production across different cell types (1). The regulation of monocyte and macrophage development to maintain homeostasis is primarily governed by M-CSF, which i...
Molekulargewicht	18.5/37.1
UniProt	P09603
Puffer	10 mM Sodium Phosphate and 50 mM Sodium Chloride, pH 8.0
Quelle	E. coli
Reinheit	97
Formulierung	Lyophilized from a sterile (0.2 micron) filtered aqueous solution

Sequenz	MEEVSEYCSH MIGSGHLQSL QRLIDSQMET SCQITFEFVD QEQLKDPVCY LKKAFLLVQD IMEDTMRFRD NTPNAIAIVQ LQELSLRLKS CFTKDYEEHD KACVRTFYET PLQLLEKVKN VFNETKNLLD KDWNIFSKNC NNSFAECSSQ GHERQSEGS
Anwendungsbeschreibung	Draw up desired volume of reconstitution buffer. Aseptically puncturing through rubber stopper with sterile needle, inject the buffer to achieve the desired concentration (0.1-0.5 mg/mL). Swirl the vial gently, do not vortex. Allow protein to rehydrate f