

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human Galectin 7 (GAL7) ELISA KBS-KBH9580

|                          |                                                                                                                                                                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human Galectin 7 (GAL7) ELISA                                                                                                                                                                                                                                  |
| Artikelnummer            | KBS-KBH9580                                                                                                                                                                                                                                                            |
| Hersteller Artikelnummer | KBH9580                                                                                                                                                                                                                                                                |
| Alternativnummer         | KBS-KBH9580-1                                                                                                                                                                                                                                                          |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                                    |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                            |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                  |
| Produktbeschreibung      | The GENLISA Human Galectin 7 (GAL7) ELISA includes features like:<br>Ready to use protocol with break-apart wells for ease of use<br>Standardisation and High Reproducibility Lot to Lot Consistency<br>Accuracy and Precision Validated against seven points for a... |
| Detektionsbereich        | 0 - 1000 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                             |
| Sensitivitaet            | 30 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                                   |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                                     |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                             |