

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Mouse Signal transducer and activator of transcription 5B (STAT5B) ELISA Kit ASC-KTE70267

|                          |                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Mouse Signal transducer and activator of transcription 5B (STAT5B) ELISA Kit                                                                 |
| Artikelnummer            | ASC-KTE70267                                                                                                                                 |
| Hersteller Artikelnummer | KTE70267                                                                                                                                     |
| Alternativnummer         | ASC-KTE70267-48,ASC-KTE70267-96                                                                                                              |
| Hersteller               | Abbkine Scientific                                                                                                                           |
| Kategorie                | Kits/Assays                                                                                                                                  |
| Applikation              | ELISA                                                                                                                                        |
| Spezies Reaktivität      | Mouse                                                                                                                                        |
| Produktbeschreibung      | This Mouse Signal transducer and activator of transcription 5B (STAT5B) ELISA Kit employs a two-site sandwich ELISA to quantitate STAT5B.... |
| Detektionsbereich        | Please inquire                                                                                                                               |
| Sensitivitaet            | Please inquire                                                                                                                               |
| Tag                      | STAT5B                                                                                                                                       |
| NCBI                     | <a href="#">20851</a>                                                                                                                        |
| UniProt                  | <a href="#">P42232</a>                                                                                                                       |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                                                                                  |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>This Mouse Signal transducer and activator of transcription 5B (STAT5B) ELISA Kit employs a two-site sandwich ELISA to quantitate STAT5B in samples. An antibody specific for STAT5B has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any STAT5B present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for STAT5B is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of STAT5B bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|