

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 488 ASC-KTL0521

Artikelname	LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 488
Artikelnummer	ASC-KTL0521
Hersteller Artikelnummer	KTL0521
Alternativnummer	ASC-KTL0521-1
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Konjugation	AbFluor 488
Produktbeschreibung	LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 488 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 488
Anwendungsbeschreibung	LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 488 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.