

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human BCAS2 (C-6His, N-T7 tag), E. coli EKL-PEH0136

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human BCAS2 (C-6His, N-T7 tag), E. coli                                                                                                                                                                                                                                |
| Artikelnummer            | EKL-PEH0136                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PEH0136                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PEH0136-10UG,EKL-PEH0136-50UG,EKL-PEH0136-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Breast Carcinoma-Amplified Sequence 2 (BCAS2) is a member of the SPF27 family. BCAS2 is a nuclear protein and widely expressed in many tissues. BCAS2 is identified as being overexpressed in various breast cancer cell lines. BCAS2 is a component of t...                       |
| Molekulargewicht         | 28.58 KDa                                                                                                                                                                                                                                                                          |
| NCBI                     | <a href="#">75934</a>                                                                                                                                                                                                                                                              |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |