

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human Carbonic Anhydrase 10 (C-6His), E. coli EKL-PEH0218

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human Carbonic Anhydrase 10 (C-6His), E. coli                                                                                                                                                                                                                          |
| Artikelnummer            | EKL-PEH0218                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PEH0218                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PEH0218-10UG,EKL-PEH0218-50UG,EKL-PEH0218-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Carbonic Anhydrase-Related Protein 10 (CA10) protein belongs to the carbonic anhydrase family of zinc metalloenzymes. It is an acatalytic member of the alpha-carbonic anhydrase subgroup. CA10 expression is detected in the adult total brain and in alm...                      |
| Molekulargewicht         | 33 KDa                                                                                                                                                                                                                                                                             |
| NCBI                     | <a href="#">9</a>                                                                                                                                                                                                                                                                  |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |