

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Goat Anti-Rabbit IgG (H&L) Highly Cross-Adsorbed - AF594 NBS-BT-CRVKZ8-100

Artikelname	Goat Anti-Rabbit IgG (H&L) Highly Cross-Adsorbed - AF594
Artikelnummer	NBS-BT-CRVKZ8-100
Hersteller Artikelnummer	BT-CRVKZ8-100
Alternativnummer	NBS-BT-CRVKZ8-100
Hersteller	Nordic BioSite
Wirt	Goat
Kategorie	Antikörper
Applikation	FC, IF
Immunogen	Rabbit IgG
Produktbeschreibung	Goat Anti-Rabbit IgG (H&L) Highly Cross-Adsorbed - AF594...
Konzentration	1 mg/ml
Puffer	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Reinheit	<p>Goat Polyclonal Secondary Antibody to Rabbit IgG (H&amp;L) have been cross-adsorbed against IgG from bovine, goat, guinea pig, hamster, horse, mouse, rat, and human. Cross-adsorption or pre-adsorption is a purification step to increase specificity of the antibody resulting in less background staining and cross-reactivity. The secondary antibody solution is passed through a column matrix containing immobilized serum proteins from potentially cross-reactive species. Only the nonspecific-binding secondary antibodies are captured in the column, and the highly specific secondaries flow through. Further passages through additional columns result in highly cross-adsorbed preparations of secondary antibody. The benefits of these extra steps are apparent in multiplexing/multicolor-staining experiments where there is potential cross-reactivity with other primary antibodies or in tissue/cell fluorescent staining experiments where there may be the presence of endogenous immunoglobulins.</p>
Antibody Type	Secondary Antibody