

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Mycobacterium ulcerans Putative O-methyltransferase MUL_4520 (MUL_4520),Escherichia Coli Preis auf Anfrage CSB-EP370746MOP

Artikelname	Recombinant Mycobacterium ulcerans Putative O-methyltransferase MUL_4520 (MUL_4520),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP370746MOP
Hersteller Artikelnummer	CSB-EP370746MOP
Alternativnummer	CSB-EP370746MOP-1, CSB-EP370746MOP-100, CSB-EP370746MOP-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A0PVW4
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MHGTDSSSDT PGQPAPSRAE LLSAHAEGSI SEDTILKTAR DRAVDIGVGA VTPAVGALLS MLAKLSGGKA VAEVGTGAGV SGLWLLSGMS DDGVLTTIDI EPEHLRVAKQ AFTEAGIGPS RTRLISGRAQ EVLTRLADDS YDLVFVDADP IDQPDYVVEG VRLLRSGGVI VVHRAALGDR AGDPAVRDAE VTAVREAARL IAEDERLTPA LVPLGDGVLA AVR D

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--