

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Xylella fastidiosa UDP-N-acetylenolpyruvoylglucosamine reductase(murB),Escherichia Coli Preis auf Anfrage CSB-EP532579XBM

Artikelname	Recombinant Xylella fastidiosa UDP-N-acetylenolpyruvoylglucosamine reductase(murB),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP532579XBM
Hersteller Artikelnummer	CSB-EP532579XBM
Alternativnummer	CSB-EP532579XBM-1, CSB-EP532579XBM-100, CSB-EP532579XBM-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0U5I2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MSRQINTPDW ILHANAPLRE LNTFHIQAQA RWLLEIIHPT ALPQALTHPH IVGLPILVLG SGSNVLFAAD PEECVLRFVN REV TILEHRI DHTLV RAGAG MAWHDLVLWS LQQGLSGLEN LALIPGTVGA CSIQNIGAYG VQVEEFVHIV EAYDQTEGKF VRLTASECEF AYRNSRFKRE PNRYLITAVE FRLPLLHELN LNYAGISEEL EALQITLPEP CDVAQAVINL RR
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.