

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Microcystis aeruginosa Glycine--tRNA ligase alpha subunit(glyQ),Escherichia Coli Preis auf Anfrage CSB-EP537114MTQ

Artikelname	Recombinant Microcystis aeruginosa Glycine--tRNA ligase alpha subunit(glyQ),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP537114MTQ
Hersteller Artikelnummer	CSB-EP537114MTQ
Alternativnummer	CSB-EP537114MTQ-1, CSB-EP537114MTQ-100, CSB-EP537114MTQ-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0JSD7
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MSLTFQAVIA KLNEFWADRG CLVAQPYDTE KGAGTMSPTH FLRAIGPEPW AVAYVEPCRR PTDGRYGENP NRFQHYYQYQ VLIKPSPDNI QEIYLDLRLV LGINPEDHDI RFVEDNWESP TLGAWGVGWE VWLDGMEITQ FTYFQQCGGI DCRPVAIEIT YGLERLAMYI QDVEAINKIY WENILYGDI FLQNEIEQCT YNFEASNPEL LFSLFSLYEQ EAKQLIDRSL VI

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--