

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Silicibacter* sp. Urease accessory protein UreF(ureF),*Escherichia Coli* Preis auf Anfrage CSB-EP624817SAAR

Artikelname	Recombinant <i>Silicibacter</i> sp. Urease accessory protein UreF(ureF), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP624817SAAR
Hersteller Artikelnummer	CSB-EP624817SAAR
Alternativnummer	CSB-EP624817SAAR-1, CSB-EP624817SAAR-100, CSB-EP624817SAAR-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q1GJQ2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	<i>E.coli</i>
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTSEATLKLM QWLSPAYPVG AFAYSHGLEG AVSEGHVRDG DTLSDWLS DL LTHGGARSDA LLLACAYRTE TAESLAEIDD TARAFPCSAE RLQETDLQGA AFCRTTAAIW DTALPSLTYP VAVGHAARLN DVALDLTLNM YLHAFSLNLV AAGQRLLSLG QTEAQRLNA CGPQLRRTVD AALSGTIDDL FGASFAADIA SMRHETQYSR IFRS

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--