

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Silicibacter* sp. ATP-dependent Clp protease ATP-binding subunit ClpX(clpX), *Escherichia Coli* Preis auf Anfrage CSB-EP625667SAAR

Artikelname	Recombinant <i>Silicibacter</i> sp. ATP-dependent Clp protease ATP-binding subunit ClpX(clpX), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP625667SAAR
Hersteller Artikelnummer	CSB-EP625667SAAR
Alternativnummer	CSB-EP625667SAAR-1, CSB-EP625667SAAR-100, CSB-EP625667SAAR-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q1GGF7
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	<i>E.coli</i>
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MATNSSGDSK NTLYCSFCGK SQHEVRKLIAGPTVFICDEC VELCMDIIRE ETKASGMKAT DGVPTPKDIC EVLDDYVIGQ ATAKRVLSVA VHNHYKRLNH AQKAGNDIEL SKSNILLIGP TGCGKTLAQ TLARILDVPF TMADATTLTE AGYVGEDVEN IILKLLQASE YNVERAQRGI VYIDEVDKIT RKSENPSITR DVSGEGVQQA LLKLMEGTVA SVPPQGGGRKH PQ

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.