

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Bacteroides fragilis N-acetylmuramic acid 6-phosphate etherase(murQ),Escherichia Coli Preis auf Anfrage CSB-EP683497BAAB

Artikelname	Recombinant Bacteroides fragilis N-acetylmuramic acid 6-phosphate etherase(murQ),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP683497BAAB
Hersteller Artikelnummer	CSB-EP683497BAAB
Alternativnummer	CSB-EP683497BAAB-1, CSB-EP683497BAAB-100, CSB-EP683497BAAB-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q5LI89
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MNSNIEKSDK PSFIKISEQP SLYDDLEKKS VREILEDINK EDQKVAIAVQ KAIPQIEKLV TQIVPRMKQG GRIFYMGAGT SGRLGVLDAS EIPPTFGMPP TLIIGLIAGG DTALRNPVEN AEDNTTRGWE ELTEHNINDK DTVIGIAASG TTPYVIGAMH AAREHGILTG CITSNPNSPM AAEDIPIEM IVGPEYVTGS SRMKSGTGQK MILNMITTSV MIQLGRVKGN KM

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.