

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Silicibacter pomeroyi* 50S ribosomal protein L3(rplC),*Escherichia Coli* Preis auf Anfrage CSB-EP683505SAAC

Artikelname	Recombinant <i>Silicibacter pomeroyi</i> 50S ribosomal protein L3(rplC), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP683505SAAC
Hersteller Artikelnummer	CSB-EP683505SAAC
Alternativnummer	CSB-EP683505SAAC-1, CSB-EP683505SAAC-100, CSB-EP683505SAAC-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q5LLU6
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MRSGIIAKKV GMTRLFMEDG KQIPVTVLQL DNLQVVAQRT TEKDGYTAVQ LGAGTPKIKR VSKAMRGHFA AQKVEPKRKL AEFRVSEDAL IEVGAEISAE HFLTGQKV DV SGTSIGKGFA GAMKRHNFGG LRASHGVSIS HRSHGSTGQC QDPGRVFKGK KMAGHMGAAR VTTQNLEVVK TDAERGLVFI KGAVPGPKSG WVTVKDAVKK KAPEGLPFFA AVKSAATEAP VA

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.