

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Silicibacter pomeroyi* Dihydrodipicolinate synthase(dapA),*Escherichia Coli* Preis auf Anfrage CSB-EP683506SAAC

Artikelname	Recombinant <i>Silicibacter pomeroyi</i> Dihydrodipicolinate synthase(dapA), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP683506SAAC
Hersteller Artikelnummer	CSB-EP683506SAAC
Alternativnummer	CSB-EP683506SAAC-1, CSB-EP683506SAAC-100, CSB-EP683506SAAC-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q5LMK7
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MFKGSMPALV TPFNGELDL EALKRLVEWQ IGEGSTGLVP VGTTGESPTL SHEEHETVVA EVVKAAAGRV PVIAGAGSNN TTEAIRFVQF AQRIGADAAL VVTPYYNRPT QRGIAHFTA LHDCAEIPV IYNIPGRSVV DMTPATMGAL AKLPRIVGVK DATGDLARVS QQRASCGADF IQLSGEDATA LGFNAHGGVG CISVTANVAP RLCAEFQQAT LAGDYAKALD YQ

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--