

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Nitratiruptor sp. Probable endonuclease 4(nfo),Escherichia Coli Preis auf Anfrage CSB-EP412359NFA

|                          |                                                                                                                                                                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Nitratiruptor sp. Probable endonuclease 4(nfo),Escherichia Coli Preis auf Anfrage                                                                                                                                                                               |
| Artikelnummer            | CSB-EP412359NFA                                                                                                                                                                                                                                                             |
| Hersteller Artikelnummer | CSB-EP412359NFA                                                                                                                                                                                                                                                             |
| Alternativnummer         | CSB-EP412359NFA-1,CSB-EP412359NFA-100,CSB-EP412359NFA-20                                                                                                                                                                                                                    |
| Hersteller               | Cusabio                                                                                                                                                                                                                                                                     |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                            |
| Spezies Reaktivität      | Bacteria                                                                                                                                                                                                                                                                    |
| Tag                      | The tag type will be determined during production process.                                                                                                                                                                                                                  |
| UniProt                  | <a href="#">A6Q1P1</a>                                                                                                                                                                                                                                                      |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                 |
| Quelle                   | E.coli                                                                                                                                                                                                                                                                      |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                 |
| Formulierung             | Lyophilized powder                                                                                                                                                                                                                                                          |
| Sequenz                  | MAKYVGAHVS ASGGVENAPK NAQAIGAKAF ALFTKNQRQW<br>KAKPLSKENI EKFNEHLECC GIEKKHILPH DSYLINLGHP EEEKRQKSLE<br>AFLDEVRRCE QLGLDKLNFH PGSHLKKISE EECLDRIADA MNETLRQSSG<br>VTLVIENTAG QGSNLGYKFE HLRYLIDKSI DKSRVGVCLD TAHMFAAGYD<br>IRTKEAYDKT MQAFDEIVGF EYLKGMHIND SKAKFASRVD RH |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|