

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Nitratiruptor sp. (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB, Escherichia Coli Preis auf Anfrage CSB-EP412375NFA

Artikelname	Recombinant Nitratiruptor sp. (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB, Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP412375NFA
Hersteller Artikelnummer	CSB-EP412375NFA
Alternativnummer	CSB-EP412375NFA-1, CSB-EP412375NFA-100, CSB-EP412375NFA-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6Q531
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MKKRLYIETL GCAMNVRDSE HIIAELTQKE DYELTQNLQE ADLILINTCS VREKPVHKL F SEIGYFNKKK KEGAKIGVCG CTASHLGEEI IKKAPYVSFV LGARNVSKIR DVIKQEKA VE VDDIDYDESTY AFSDFRTSPY KAFINISIGC DKKCTFCIVP NTRGEEISIP PELILQEVRR AVDTGAKEIF LLGQNVNNYG RSFSDKARKV DFTDLLRMVS EIDGVRRIRF TS

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--