

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Sulfurovum sp. Imidazoleglycerol-phosphate dehydratase(hisB),Escherichia Coli Preis auf Anfrage CSB-EP412390SUZ

Artikelname	Recombinant Sulfurovum sp. Imidazoleglycerol-phosphate dehydratase(hisB),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP412390SUZ
Hersteller Artikelnummer	CSB-EP412390SUZ
Alternativnummer	CSB-EP412390SUZ-1,CSB-EP412390SUZ-100,CSB-EP412390SUZ-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6Q9M2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MIEVTRETKE TQISVKLNLY GKGQAKIDTG VGFFDHMLEA FTKHAHIDME VTCTGDTHID DHHTVEDVGI VIGQALREAI YPVKSIERFG NATVVMDEAS VTCIDIDLSNR GFLVFELPIG GKVGEFDVEL VEEFFRAFAF NLPLTLHLIY NRGKNKHII EAAFKALAVS LRRAVAINEN AGIPSTKGV

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--