

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Klebsiella pneumoniae subsp. pneumoniae Autoinducer 2 import ATP-binding protein LsrA(IsrA),Escherichia Coli Preis auf Anfrage CSB-EP412572KAX

Artikelname	Recombinant Klebsiella pneumoniae subsp. pneumoniae Autoinducer 2 import ATP-binding protein LsrA(IsrA),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP412572KAX
Hersteller Artikelnummer	CSB-EP412572KAX
Alternativnummer	CSB-EP412572KAX-1,CSB-EP412572KAX-100,CSB-EP412572KAX-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6TEB8
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MMKPLLEARQ ICKQFSGVAV LKGIDFTLGA GQVHALMGGN GAGKSTLMKI IAGVETPD SG ELTIGERAFA RLNPALAHQL GIYLV PQEPM LFPNLSVREN ILFRLPKRAD TAARLQDKLQ QLNCQINLDA SASTLEVADQ QMVEILRGLM REARILILDE PTASLTPGET ERLFSQIRAL QALDVGIVFI SHKLPEIRQL ASHISVMRDG AVVLSGETAA YRDEQLINAM TP

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--