

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Dinoroseobacter shibae (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB, Escherichia Coli Preis auf Anfrage CSB-EP421812DIB

Artikelname	Recombinant Dinoroseobacter shibae (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB, Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP421812DIB
Hersteller Artikelnummer	CSB-EP421812DIB
Alternativnummer	CSB-EP421812DIB-1, CSB-EP421812DIB-100, CSB-EP421812DIB-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A8LSE7
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MADPKKLYIK TYGCQMNVDYD SERMAEAMGG EGYVQTDRAE DADMILLNTC HIREKAAEKV YSELGRFKPL KAEKPD LKIG VAGCVAQAEG AEIMRRQPLV DLVVG PQAYH RLP ELTARAA TGAKALD TDF PEEDKFDHLA ARPKAKRGPT AFLTVQEGCD KFC AFCVVPY TRGA EVSRPA ARVLTEARDL VERGVREITL LGQNVNAYHG HARGLAGLIW DLAEIDGLER IR

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--