

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Nitrosopumilus maritimus 50S ribosomal protein L2P(rpl2p),Escherichia Coli Preis auf Anfrage CSB-EP427857NFB

|                          |                                                                                                                                                                                                                                                                               |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Nitrosopumilus maritimus 50S ribosomal protein L2P(rpl2p),Escherichia Coli Preis auf Anfrage                                                                                                                                                                      |
| Artikelnummer            | CSB-EP427857NFB                                                                                                                                                                                                                                                               |
| Hersteller Artikelnummer | CSB-EP427857NFB                                                                                                                                                                                                                                                               |
| Alternativnummer         | CSB-EP427857NFB-1, CSB-EP427857NFB-100, CSB-EP427857NFB-20                                                                                                                                                                                                                    |
| Hersteller               | Cusabio                                                                                                                                                                                                                                                                       |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                              |
| Spezies Reaktivität      | Bacteria                                                                                                                                                                                                                                                                      |
| Tag                      | The tag type will be determined during production process.                                                                                                                                                                                                                    |
| UniProt                  | <a href="#">A9A141</a>                                                                                                                                                                                                                                                        |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                   |
| Quelle                   | E.coli                                                                                                                                                                                                                                                                        |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                   |
| Formulierung             | Lyophilized powder                                                                                                                                                                                                                                                            |
| Sequenz                  | MGKRPLVRRR GRGGNQFRST STGKVGTKAN YPRFPLSEQH EGEIIDLVHE<br>RGREAPLAKV RFEDGSVSFV PAVLGAKVGE TLQFGLKSKI EKGNVISIQN<br>IPDGTIVCNI EKHF GDGGAI VKSAGTDATV FSHGDEGVTV KLPSGKFTTL<br>NPKNRAMIGT LAGGGASERP FMSAGGKWRN FKAKGKKYPI<br>VRGVAQAAYV HPHGGGRHQH VGQSSTVSRD TPPGAKVGSII AA |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|