

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Burkholderia multivorans (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB, Escherichia Coli Preis auf Anfrage CSB-EP432560BPS

Artikelname	Recombinant Burkholderia multivorans (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB, Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP432560BPS
Hersteller Artikelnummer	CSB-EP432560BPS
Alternativnummer	CSB-EP432560BPS-1, CSB-EP432560BPS-100, CSB-EP432560BPS-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A9AFF8
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTKKVYVKTF GCQMNEYDSD KMVDVLNAAE GLEKTDSPED ADIILFNTCS VREKAQEKVF SDLGRVRELK EAKPDLLIGV GGCVASQEGA SIVSRAPYVD LVFGPQTLHR LPQMIDARRT SGRAQVDITF PEIEKFDHLP PARVEGPTAF VSIMEGCSKY CSYCVVPYTR GDEVSRPLDD VLTEVAGLAD QGVREVTLLG QNVNAYRGAL TAGSSEVADF ATLIEYVADI PG

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--