

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Burkholderia phymatum Formamidopyrimidine-DNA glycosylase(mutM),Escherichia Coli Preis auf Anfrage CSB-EP450564BXT

Artikelname	Recombinant Burkholderia phymatum Formamidopyrimidine-DNA glycosylase(mutM),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP450564BXT
Hersteller Artikelnummer	CSB-EP450564BXT
Alternativnummer	CSB-EP450564BXT-1, CSB-EP450564BXT-100, CSB-EP450564BXT-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	B2JCP3
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	PELPEVEVT RRGIEPYVAG RRVERVDVRT PALRWPIPPG FARLLHGRLV RKVERRGKYL LFEIDEGWFI VHLGMTGTLR VLRNVPHPPA AAKHDHVDWI FDDFILRFRD PRRFGAVLWH PRSAGDVLDH PLLADLGVEP FAPSFSGALL HRKTRGRKVS VKQALLAGEI VVGVGNIYAS ESLFRAGIRP TTPAGRISLV RYDLLADAVR VTLAAAIEKG GSTLRDFVGS NGE

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--