

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Danio rerio RING finger protein 141(rnf141),Escherichia Coli Preis auf Anfrage CSB-EP019839DIL

Artikelname	Recombinant Danio rerio RING finger protein 141(rnf141),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP019839DIL
Hersteller Artikelnummer	CSB-EP019839DIL
Alternativnummer	CSB-EP019839DIL-1, CSB-EP019839DIL-100, CSB-EP019839DIL-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Fish
Tag	The tag type will be determined during production process.
UniProt	Q6IV56
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MGQQLSGQAV TRLPEKLIKH VGLVRDSGYL TYDEFLGRVA ELNDVTAKLA SGQKKHLLFE VQPGSDSSAL WKVAVRIVCT KINKEDGMVE ASRIMNLYQF IQLYKDITSQ AAEVFSSNVA AEGSSEDTCQ ASMWMGRVKQ LTDEEECCIC MDGKADLILP CAHSFCQKCI DKWSGQSRNC PVCRIQVTAA NESWVMSDAP TGEDVAGYIL NLADEAGHPH RP

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--