

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Macaca fascicularis RING finger protein 32(RNF32), Escherichia Coli Preis auf Anfrage CSB-EP019885MOV

Artikelname	Recombinant Macaca fascicularis RING finger protein 32(RNF32), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP019885MOV
Hersteller Artikelnummer	CSB-EP019885MOV
Alternativnummer	CSB-EP019885MOV-1, CSB-EP019885MOV-100, CSB-EP019885MOV-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Monkey
Tag	The tag type will be determined during production process.
UniProt	Q4R5T4
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MLKNKGHSSK KDNLAVNAVA LQDHILRDLQ LRNLSVADHS KTQVQKKENK SLKRDTKAII DTGLKKTMCQ PKLEDSEKEY VLDPKPPPLT LAQKLGLIGP PPPQLSSDEW EKVQQRSLQ GDSVQPCPIC KEEFELRPQV LLSCSHVFHR ACLQAF EKFT NKKTCPLCRK NQYQTRVIHG GARLFRIKCV TRIQAYWRGY VVRKWYRNLR ETVPPTDAKL RKKFFEKKFT EI

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--