

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Nematostella vectensis Ubiquitin-conjugating enzyme E2 S (v1g237158), Escherichia Coli Preis auf Anfrage CSB-EP025477NGO

Artikelname	Recombinant Nematostella vectensis Ubiquitin-conjugating enzyme E2 S (v1g237158), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP025477NGO
Hersteller Artikelnummer	CSB-EP025477NGO
Alternativnummer	CSB-EP025477NGO-1, CSB-EP025477NGO-100, CSB-EP025477NGO-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A7SE05
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MATTSSNVEN FSPQIIKQVA REIHGLTNDP PEGIKVFSND EDITDIQASI EGPTGTPYEG GIFKIKLVLG KDFPAAPPKG FFLTKIFHPN VAKNGEICVN TLKKDWKPD L GIKQVLLTVK CLLIVPNPES ALNEEAGKLL LERYDDYSKR AKMFTEIHAK LSASSNNIS EGQQESLP GK KRVAVNEKMC DKKKKDKKRA LKRL

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.