

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

**Recombinant Salmonella typhimurium
Phosphoribosylformylglycinamide cyclo-ligase(purM), Escherichia Coli
Preis auf Anfrage
CSB-EP301886SXB**

Artikelname	Recombinant Salmonella typhimurium Phosphoribosylformylglycinamide cyclo-ligase(purM), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP301886SXB
Hersteller Artikelnummer	CSB-EP301886SXB
Alternativnummer	CSB-EP301886SXB-1, CSB-EP301886SXB-100, CSB-EP301886SXB-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	P74883
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	TDKTSLSYK DAGVDIDAGN ALVDRIKGVV KKTRRPEVMG GLGGFGALCA LPQKYREPV L VSGTDGVGTK LRLAMD LKRH DAIGIDL VAM CVNDLVVQGA EPLFFLDYYA TGKLDVDTAA SVINGIAEGC LQSGCALVGG ETAEMPGMYH GEDYDVAGFC VGVVEKSEII DGSRVAEGDV LIALGSSGPH SNGYSLVRKI IDVSGCDPQT TLLEGKPLAD HLLPTRIYV KSV
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.