

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Tobacco mosaic virus Capsid protein (CP) CSB-EP304020TFF

Artikelname	Recombinant Tobacco mosaic virus Capsid protein (CP)
Artikelnummer	CSB-EP304020TFF
Hersteller Artikelnummer	CSB-EP304020TFF
Alternativnummer	CSB-EP304020TFF-1,CSB-EP304020TFF-100,CSB-EP304020TFF-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Molekulargewicht	24.4 kDa
Tag	C-terminal 6xHis-tagged
UniProt	P69687
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	E.coli
Expression System	2-159aa
Reinheit	Greater than 95% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder
Sequenz	SYSITTPSQFVFLSSAWADPIELINLCTNALGNQFQTQQARTVVQRQFSEVWK PSPQVTVRFPDSDFKVYRYNAVLDPDPLVTALLGAFDTRNRIIEVENQANPTTAET LDATRRVDDATVAIRSAINNLLIVELIRGTGSYNRSSFESSSGLVWTSGPAT

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.