

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Hirudo medicinalis* Metallo-carboxypeptidase inhibitor, *Escherichia Coli* Preis auf Anfrage CSB-EP305666HSM

Artikelname	Recombinant <i>Hirudo medicinalis</i> Metallo-carboxypeptidase inhibitor, <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP305666HSM
Hersteller Artikelnummer	CSB-EP305666HSM
Alternativnummer	CSB-EP305666HSM-1, CSB-EP305666HSM-100, CSB-EP305666HSM-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Worm
Tag	The tag type will be determined during production process.
UniProt	P81511
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	SHTPD ESFLCYQPDQ VCCFICRGAA PLPSEGECPNP HPTAPWCREG AVEWVPYSTG QCRTTCIPYV E

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--